


# STOOL WITH BACKREST



A4003.1025B



A3202.1025B



A3069.1025B



A2001.1025B



A5069.1025B

## General Specifications

### DESCRIPTION

Elegant stool with a backrest ergonomically designed with a flat round seat.

Equipped with 3 levers to adjust height, backrest and seat inclination.

### UPHOLSTERY

Spradling PU

UPHOLSTERY SIZE  
Ø37 x 7,5 cm

BASE DIAMETER  
64 cm

UP AND DOWN (SEAT)  
51 / 71 cm

### PACKAGING SIZE

59 x 30 x 60 cm

N.W. KG  
8

G.W. KG  
10

### OPTIONS

Lower gas lift:  
Ref. 30004

UP AND DOWN  
47/67 cm