

LATIS

2-section table with a very stable metallic structure. Equipped with 1 motor that controls the table height. Backrest adjustable by gas lift. Removable and adjustable headrest. 2-wheel system to facilitate its mobility. White PU upholstery of high quality that eases cleaning.

Camilla de 2 cuerpos con estructura metálica muy estable. Equipada con 1 motor que permite regular la altura. Respaldo regulable mediante pistón de gas. Reposacabezas regulable y extraíble. Incluye 2 ruedas que facilitan su desplazamiento. Tapizado blanco en PU de alta calidad y fácil de limpiar.

Table à 2 plans avec une structure métallique très stable. Équipée d'1 moteur qui règle la hauteur. Dossier réglable par vérin à gaz. Tête réglable et amovible. Système à 2 roulettes pour faciliter son déplacement. Revêtement en PU blanc de haute qualité et facile à nettoyer.

Sehr stabile zweiteilige Liege. Sie ist mit einem Motor ausgestattet, um die Höhe zu regulieren. Mittels einer Gasdruckfeder ist auch die Neigung der Rückenlehne regulierbar. Ausziehbare und verstellbare Kopfstütze. Struktur mit 2 Rädern, um die Verschiebung zu vereinfachen. Hochwertige weiße PU-Polsterung, leicht abwaschbar.

Letino a 2 snodi molto stabile. Dotato di 1 motore che regola l'altezza. L'inclinazione dello schienale è regolabile con un pistone a gas. Poggiatesta estraibile con foro facciale. Struttura con 2 ruote per un facile spostamento. Rivestimento in PU bianco di grande qualità e facile da pulire.

Electric Massage Bed Ref. 2210A.1.A26

Product size: 204x68x68 cm
Weight: 50 Kg
Upholstery: PU
Motors: 1

- 1 Angle/Extension/Breather hole
- 2 Angle 0° ~ 55°
- 3 Up and Down 68/92 cm
- 4 Lock/unlock wheels

OPTIONS

- Bed cover
- 5 Hand controller Ref. HC001
- 6 Foot controller Ref. FC001
- Paper roll holder

PACKING SIZE cm

180x72x49

N.W. kg

50

G.W. kg

58

MAX WEIGHT

150 kg



LATIS

